





LENS DESIGN

PORTFOLIO

*"Personalized, digital free-form lenses
made to every patient's specific needs,
using the latest innovations in lens
optimization technology"*

IOT Digital Ray-Path 2



Progressive Lenses

Endless Progressive

BEST

"Optimized to ensure impeccable vision quality"

Fully personalized progressive lens with advanced design that considers the wearer's accommodative ability, ensuring exceptional clarity and comfort at all distances.

- Precise, comfortable focus for all working distances
- Near elimination of peripheral blur
- Superior visual quality with digital devices
- Higher image stability, reduced swim effect
- Improved peripheral acuity in distance zone

MFH: 14,15,16,17,18 mm

Steady Progressive

BETTER

"High quality and affordable progressive lens"

Non-compensated free-form progressive lens offering a larger field of vision and greater comfort than others in its class.

- Higher image stability, reduced swim effect
- Improved peripheral acuity in distance zone
- Good balance between near and distance fields
- Variable inset for improved binocular vision

MFH: 14,15,16,17,18 mm

Essential Progressive

GOOD

"Entry-level digital progressive lens"

Non-compensated progressive lens made with free-form technology, offering superior vision to a conventional lens.

- Good balance between near and distance fields
- Variable inset for improved binocular vision and wider near area

MFH: 14,15,16,17,18 mm

Classic Adapt

BASIC

"Basic conventional progressive lens"

Conventional progressive lens offering easy and comfortable adaptation.

- Balanced near and distance fields
- Fixed inset

MFH: 14, 17 mm

Specialty lenses · Single Vision Lenses

Endless Office

OCCUPATIONAL

Personalized occupational lens with maximum near and intermediate vision zones for superb vision in the workspace.

WORKSTATION

35 cm – 2 m / 14 in – 6.5 ft

Perfect for wearers needing excellent desk-level vision while engaging with people at conversational distance.

- Improved postural ergonomics, avoiding unnecessary head movements
- Comfortable and precise focusing with electronic devices
- Near elimination of peripheral blur
- Precise focus for all working distances in any direction of gaze

MFH: 14, 18 mm

ROOM

35 cm – 4 m / 14 in – 13.1 ft

Perfect for wearers moving around a larger workspace who benefit from additional distance viewing.

- Excellent dynamic vision, easy near-to-intermediate transition
- Near elimination of peripheral blur

MFH: 14, 18 mm

SCREEN

35 cm – 1.3 m / 14 in – 4.2 ft

Perfect for wearers spending most of their time focusing on nearby objects, especially a computer screen.

- More relaxed vision with less accommodative effort
- Near elimination of peripheral blur

MFH: 16, 18 mm

Endless

Single Vision

PREMIUM

Premium personalized single-vision lens designed for everyday use, with impeccable visual quality, clarity, and comfort.

- Consistently excellent visual quality, including high prescriptions and wrapped frames
- Near elimination of peripheral blur
- Excellent dynamic vision with minimal distortion
- Superior visual quality with digital devices

Endless

Anti-Fatigue SV

FATIGUE RELIEF

Personalized single-vision lens featuring a unique power boost to combat visual fatigue, featuring 3 available power boosts.

Power Boosts Available:

0.50 D · 0.75 D · 1.00 D

- More relaxed vision with less accommodative effort
- Excellent distance & peripheral vision

GROUNDBREAKING TECHNOLOGY

IOT Digital Lens Optimization · Core Technology Principles

99.5%

of gaze directions now optimized

Spatial Vision

Optimizes the eye's ability to focus on near and far objects with unparalleled sharpness perception and extended focus areas. Enables smooth and rapid transitions across visual zones, especially when looking down to focus on nearby objects.

WITHOUT SPATIAL VISION



IOT Digital Ray-Path 2 · US Patent 11,609,437 B2

WITH SPATIAL VISION



IOT Digital Ray-Path 2

Maximizes the clear viewing area and significantly reduces perceived peripheral blur, achieving an absolute reduction in defocus and increasing visual performance of the lens.

Eye Focus Profile

Point-by-point optimization over the entire lens surface provides precise vision at every distance and in every direction of gaze, for a more natural visual experience.

Ray Tracing

Calculates the back surface using a pure geometrical method, producing lenses with flexible designs and variable corridor lengths and insets.

Binocular Vision

Significantly reduces swim effect when moving. Virtually eliminates image distortion and enhances image stability for a comfortable and clear visual experience.

Visual Stability

Virtually eliminates image distortion and enhances image stability for a comfortable and clear visual experience when the wearer is in motion.

Digital Surfacing

Enables flexible designs and variable corridor lengths through digital surfacing methods, allowing highly personalized lens geometry.

Digital Vision

Designed to alleviate eyestrain caused by frequent digital device use, providing optimal visual support for screen-based work.

These principles underpin IOT's cutting-edge lens optimization: Camber Technology · IOT Digital Ray-Path 2 · Steady Plus Methodology